

IFN Gamma / Interferon Gamma Mouse anti-Bovine Monoclonal (RPE) (CC302) Antibody - LS-C188184 - LSBio

CatalogID:	LS-C188184
Target:	interferon, gamma (IFNG)
Synonyms:	IFNG Antibody, IFG Antibody, Immune interferon Antibody, IFI Antibody, Interferon gamma Antibody, IFN-gamma Antibody, Interferon, gamma Antibody
Family / Subfamily:	Interferon / not assigned-Interferon
Host	IFNG antibody was produced in Mouse
Clonality:	Monoclonal
Isotype:	IgG1
Clone Name:	CC302
Conjugations:	R. Phycoerythrin (RPE)
Immunogen Species:	IFN Gamma / Interferon Gamma antibody was raised against Bovine
Reactivity:	Bovine, Human, Dog, Dolphin, Ferret, Goat, Horse, Mink, Pig, Sheep, Fin whale
Predicted Reactivity:	Mustelid
Purification:	Protein G purified
Reconstitution:	Distilled water
Presentation:	PBS, 0.09% sodium azide, 1% BSA, 5% sucrose
Recommended Storage:	Long term: -20°C; Short term: +4°C; Avoid freeze-thaw cycles. Protect from light.
Uses:	Flow Cytometry (1:1 - 1:10) (Optimal dilution to be determined by the researcher)
Size:	100 tst
Requested From:	Japan

Laboratory Reagent For In Vitro Research Use Only

Not for resale without prior written consent from LifeSpan BioSciences, Inc.

Created on 9/26/2014

© 2014 LifeSpan BioSciences