

IL4 Mouse anti-Bovine Monoclonal (RPE) (CC303) Antibody - LS-C188186 - LSBio

CatalogID:	LS-C188186
Target:	interleukin 4 (IL4)
Synonyms:	IL4 Antibody, B cell growth factor 1 Antibody, B-cell stimulatory factor 1 Antibody, B_cell stimulatory factor 1 Antibody, BCGF1 Antibody, BSF-1 Antibody, Binetrakin Antibody, BCGF-1 Antibody, Interleukin 4 Antibody, Interleukin-4 Antibody, Pitrakinra Antibody, BSF1 Antibody, IL-4 Antibody
Family / Subfamily:	Hormone / not assigned-Hormone
Host	IL4 antibody was produced in Mouse
Clonality:	Monoclonal
Isotype:	IgG2a
Clone Name:	CC303
Conjugations:	R. Phycoerythrin (RPE)
Immunogen Species:	IL4 antibody was raised against Bovine
Reactivity:	Bovine, Dog, Dolphin, Goat, Mink, Pig, Sheep, Fin whale, Mustelid
Purification:	Protein G purified
Reconstitution:	Distilled water
Presentation:	PBS, 0.09% sodium azide, 1% BSA, 5% sucrose
Recommended Storage:	Long term: -20°C; Short term: +4°C; Avoid freeze-thaw cycles. Protect from light.
Uses:	Flow Cytometry (1:1 - 1:5) (Optimal dilution to be determined by the researcher)
Size:	100 tst
Requested From:	Japan

Laboratory Reagent For In Vitro Research Use Only

Not for resale without prior written consent from LifeSpan BioSciences, Inc.

Created on 9/26/2014

© 2014 LifeSpan BioSciences