

CD45RB Rat anti-Mouse Monoclonal (G1/G1-4) Antibody - LS-C189683 - LSBio

|                             |                                                                                                                                                 |
|-----------------------------|-------------------------------------------------------------------------------------------------------------------------------------------------|
| <b>CatalogID:</b>           | LS-C189683                                                                                                                                      |
| <b>Target:</b>              | CD45RB                                                                                                                                          |
| <b>Host</b>                 | CD45RB antibody was produced in Rat                                                                                                             |
| <b>Clonality:</b>           | Monoclonal                                                                                                                                      |
| <b>Isotype:</b>             | IgG2a,k                                                                                                                                         |
| <b>Clone Name:</b>          | G1/G1-4                                                                                                                                         |
| <b>Immunogen Species:</b>   | CD45RB antibody was raised against Mouse                                                                                                        |
| <b>Immunogen:</b>           | CD45RB antibody was raised against cloned mouse Th2 cell lines.                                                                                 |
| <b>Specificity:</b>         | Mouse CD45RB                                                                                                                                    |
| <b>Reactivity:</b>          | Mouse                                                                                                                                           |
| <b>Purification:</b>        | Immunoaffinity purified                                                                                                                         |
| <b>Presentation:</b>        | PBS, pH 7.2, 0.09% sodium azide                                                                                                                 |
| <b>Recommended Storage:</b> | May be stored at 4°C for short-term only. Aliquot to avoid repeated freezing and thawing. Store at -20°C. Aliquots are stable for at 12 months. |
| <b>Uses:</b>                | IHC, Immunoprecipitation, Flow Cytometry (Optimal dilution to be determined by the researcher)                                                  |
| <b>Size:</b>                | 500 µg                                                                                                                                          |
| <b>Concentration:</b>       | 0.5 mg/ml                                                                                                                                       |
| <b>Requested From:</b>      | Japan                                                                                                                                           |

Laboratory Reagent For In Vitro Research Use Only

Not for resale without prior written consent from LifeSpan BioSciences, Inc.

Created on 9/27/2014

© 2014 LifeSpan BioSciences