

CD1A Mouse anti-Human Monoclonal (aa1-300) (PE) (1D4B) Antibody - LS-C189779 - LSBio	
CatalogID:	LS-C189779
Target:	CD1a molecule
Synonyms:	CD1A Antibody, CD1 Antibody, CD1A antigen, a polypeptide Antibody, CD1a molecule Antibody, FCB6 Antibody, HTA1 Antibody, HTa1 thymocyte antigen Antibody, R4 Antibody, T6 Antibody, CD1a antigen Antibody, Cluster of differentiation 1 A Antibody
Host	CD1A antibody was produced in Mouse
Clonality:	Monoclonal
Isotype:	IgG1
Clone Name:	1D4B
Conjugations:	Phycoerythrin (PE)
Immunogen Species:	CD1A antibody was raised against Human
Antigen Type:	Recombinant protein
Immunogen:	CD1A antibody was raised against recombinant corresponding to aa1-300 from human CD1a expressed in NSO cell line.
Specificity:	Human CD1A.
Epitope:	aa1-300
Reactivity:	Human
Purification:	Affinity purified
Presentation:	Saline, 0.5% BSA, 0.1% sodium azide
Recommended Storage:	Maybe stored at 4°C before opening. DO NOT FREEZE! Stable at 4°C as an undiluted liquid. Stable for 12 months at 4°C. Sensitive to light
Uses:	Flow Cytometry (Optimal dilution to be determined by the researcher)
Size:	100 tst
Requested From:	Japan
Laboratory Reagent For In Vitro Research Use Only	
Not for resale without prior written consent from LifeSpan BioSciences, Inc.	
Created on 9/27/2014	
© 2014 LifeSpan BioSciences	