

LY6G / Gr-1 Rat anti-Mouse Monoclonal (PE) (GF1) Antibody - LS-C189923 - LSBio	
CatalogID:	LS-C189923
Target:	lymphocyte antigen 6 complex, locus G (LY6G)
Synonyms:	LY6G Antibody, Ly-6G Antibody, Gr-1 Antibody, Gr1 Antibody, Ly6g Antibody
Host	LY6G antibody was produced in Rat
Clonality:	Monoclonal
Isotype:	IgG2b
Clone Name:	GF1
Conjugations:	Phycoerythrin (PE)
Immunogen Species:	LY6G / Gr-1 antibody was raised against Mouse
Antigen Type:	Cells
Immunogen:	LY6G / Gr-1 antibody was raised against purified BALB/c mouse neutrophils.
Specificity:	Mouse Ly-6G and Ly-6C
Reactivity:	Mouse
Purification:	Protein G purified
Presentation:	PBS, 1% BSA, 0.02% sodium azide
Recommended Storage:	Maybe stored at 4°C before opening. DO NOT FREEZE! Stable at 4°C as an undiluted liquid. Stable for 12 months at 4°C. Sensitive to light
Usage Summary:	Functional Studies: Neutrophil depletion. Mice were treated with NIMP-R1.
Uses:	IHC - Paraffin, IHC - Frozen, Flow Cytometry, ELISA, Blocking, Functional Assay (Optimal dilution to be determined by the researcher)
Size:	25 tst
Requested From:	Japan
Laboratory Reagent For In Vitro Research Use Only	
Not for resale without prior written consent from LifeSpan BioSciences, Inc.	
Created on 9/27/2014 © 2014 LifeSpan BioSciences	