

FITC Mouse Monoclonal (FITC) (NI 239) Antibody - LS-C191526 - LSBio	
CatalogID:	LS-C191526
Target:	Fluorescein Isothiocyanate (FITC)
Host	Fluorescein Isothiocyanate (FITC) antibody was produced in Mouse
Clonality:	Monoclonal
Isotype:	IgG1,k
Clone Name:	NI 239
Conjugations:	Fluorescein (FITC)
Antigen Type:	Purified protein
Immunogen:	FITC antibody was raised against highly purified fluorescein isothiocyanate isomer 1.
Specificity:	Fluorescein isothiocyanate.
Purification:	Ascites
Reconstitution:	Sterile distilled water. Possible additional volumes for resuspension: 500 µl
Presentation:	Lyophilized
Recommended Storage:	Lyophilized powder may be stored at -20°C. Stable for 12 months at -20°C. Aliquot to avoid repeated freezing and thawing. Store at -20°C. Reconstituted product is stable for 12 months at -20°C.
Uses:	Western blot, ELISA (Optimal dilution to be determined by the researcher)
Size:	500 µg
Requested From:	Japan
Laboratory Reagent For In Vitro Research Use Only Not for resale without prior written consent from LifeSpan BioSciences, Inc. Created on 9/27/2014 © 2014 LifeSpan BioSciences	