

F8 / FVIII / Factor VIII Sheep anti-Human Polyclonal (FITC) Antibody - LS-C193839 - LSBio	
CatalogID:	LS-C193839
Target:	coagulation factor VIII, procoagulant component (F8)
Synonyms:	F8 Antibody, AHF Antibody, Antihemophilic factor Antibody, Coagulation factor VIIIc Antibody, Coagulation factor VIII Antibody, DXS1253E Antibody, Factor VIII F8B Antibody, FVIII Antibody, HEMA Antibody, F8B Antibody, Factor VIII Antibody, Procoagulant component Antibody, F8C Antibody
Host	F8 antibody was produced in Sheep
Clonality:	Polyclonal
Isotype:	IgG
Conjugations:	Fluorescein (FITC)
Immunogen Species:	F8 / FVIII / Factor VIII antibody was raised against Human
Immunogen:	F8 / FVIII / Factor VIII antibody was raised against human Factor VIII purified from concentrate.
Specificity:	Human Factor VIII
Reactivity:	Human
Purification:	Affinity purified
Presentation:	PBS, 0.1% BSA, 0.1% sodium azide, pH 7.4
Recommended Storage:	Store at 4°C. Do not freeze. Protect from light.
Uses:	Undetermined
Size:	100 µg
Requested From:	Japan
Laboratory Reagent For In Vitro Research Use Only	
Not for resale without prior written consent from LifeSpan BioSciences, Inc.	
Created on 9/27/2014 © 2014 LifeSpan BioSciences	