

Luciferase Goat anti-Firefly Polyclonal (FITC) Antibody - LS-C195078 - LSBio	
<b>CatalogID:</b>	LS-C195078
<b>Target:</b>	Luciferase
<b>Host</b>	Luciferase antibody was produced in Goat
<b>Clonality:</b>	Polyclonal
<b>Conjugations:</b>	Fluorescein (FITC)
<b>Immunogen Species:</b>	Luciferase antibody was raised against Firefly
<b>Immunogen:</b>	Luciferase antibody was raised against highly purified firefly luciferase.
<b>Specificity:</b>	Firefly Luciferase
<b>Reactivity:</b>	Firefly
<b>Purification:</b>	Affinity purified
<b>Presentation:</b>	PBS, pH 7.4, 0.2% BSA, 0.05% sodium azide
<b>Recommended Storage:</b>	Short term 4°C, long term aliquot and store at -20°C, avoid freeze thaw cycles.
<b>Uses:</b>	ICC (1:10), ELISA (1:10000) (Optimal dilution to be determined by the researcher)
<b>Size:</b>	500 µl
<b>Concentration:</b>	50 µg/ml
<b>Requested From:</b>	Japan
<b>Laboratory Reagent For In Vitro Research Use Only</b>	
Not for resale without prior written consent from LifeSpan BioSciences, Inc.	
Created on 9/27/2014	
© 2014 LifeSpan BioSciences	