

Hemoglobin Mouse anti-Human Monoclonal (PE) (70050) Antibody - LS-C196104 - LSBio

CatalogID:	LS-C196104
Target:	Hemoglobin
Host	Hemoglobin antibody was produced in Mouse
Clonality:	Monoclonal
Isotype:	IgG1,k
Clone Name:	70050
Conjugations:	Phycoerythrin (PE)
Immunogen Species:	Hemoglobin antibody was raised against Human
Antigen Type:	Synthetic peptide
Immunogen:	Hemoglobin antibody was raised against a synthetic peptide sequence unique to the gamma chain of embryonic human Hemoglobin.
Specificity:	Hemoglobin gamma. Hb alpha - chain < 0.1%
Reactivity:	Human
Purification:	Purified IgG
Presentation:	PBS, pH 7.4, 1.0% BSA, 0.01% sodium azide
Recommended Storage:	Store at 4°C. Do not freeze.
Uses:	Undetermined
Size:	100 µg
Requested From:	Japan

Laboratory Reagent For In Vitro Research Use Only

Not for resale without prior written consent from LifeSpan BioSciences, Inc.

Created on 9/27/2014

© 2014 LifeSpan BioSciences