

PECAM-1 / CD31 Mouse anti-Human Monoclonal (FITC) (HC1/6) Antibody - LS-C196260 - LSBio	
CatalogID:	LS-C196260
Target:	platelet/endothelial cell adhesion molecule 1 (PECAM1)
Synonyms:	PECAM1 Antibody, CD31 Antibody, CD31 antigen Antibody, EndoCAM Antibody, PECA1 Antibody, PECAM-1 Antibody, GPIIA Antibody
Host	PECAM1 antibody was produced in Mouse
Clonality:	Monoclonal
Isotype:	IgG1
Clone Name:	HC1/6
Conjugations:	Fluorescein (FITC)
Immunogen Species:	PECAM-1 / CD31 antibody was raised against Human
Antigen Type:	Cells
Immunogen:	PECAM-1 / CD31 antibody was raised against pMA-treated U937 cells.
Specificity:	Human PECAM / CD31
Reactivity:	Human
Purification:	Affinity purified
Presentation:	PBS, 10 mM sodium azide, 0.1% BSA
Recommended Storage:	Store at 4°C. Do not freeze. Protect from light.
Uses:	Flow Cytometry (Optimal dilution to be determined by the researcher)
Size:	100 tst
Requested From:	Japan
Laboratory Reagent For In Vitro Research Use Only	
Not for resale without prior written consent from LifeSpan BioSciences, Inc.	
Created on 9/27/2014	
© 2014 LifeSpan BioSciences	