


OVOL1 Mouse anti-Human Monoclonal (1B3) Antibody - LS-C197830 - LSBio

CatalogID:	LS-C197830
Target:	ovo-like zinc finger 1 (OVOL1)
Synonyms:	OVOL1 Antibody, HOVO1 Antibody, Ovo-like 1(Drosophila) Antibody, Ovo homolog-like 1 Antibody
Host	OVOL1 antibody was produced in Mouse
Clonality:	Monoclonal
Isotype:	IgG1,k
Clone Name:	1B3
Immunogen Species:	OVOL1 antibody was raised against Human
Antigen Type:	GST fusion protein
Immunogen:	OVOL1 antibody was raised against oVOL1 (AAH59408, aa1-267 full-length recombinant protein with GST tag. MW of the GST tag alone is 26 KDa.
Specificity:	Human OVOL1
Reactivity:	Human
Purification:	Protein A purified
Presentation:	PBS, pH 7.2
Recommended Storage:	Store at -20°C or lower. Aliquot to avoid repeated freezing and thawing.
Usage Summary:	Sandwich ELISA: Recombinant protein.
Uses:	ELISA (Optimal dilution to be determined by the researcher)
Size:	100 µg

ELISA Image:



Detection limit for recombinant GST tagged OVOL1 is 3 ng/ml as a capture antibody.

Requested From: Japan

Laboratory Reagent For In Vitro Research Use Only

Not for resale without prior written consent from LifeSpan BioSciences, Inc.

Created on 9/27/2014

© 2014 LifeSpan BioSciences