

Mitocryptide-1 Rabbit anti-Human Polyclonal (Cy5) Antibody - LS-C203528 - LSBio

CatalogID:	LS-C203528
Target:	Microcryptide-1
Host	Microcryptide-1 antibody was produced in Rabbit
Clonality:	Polyclonal
Isotype:	IgG
Conjugations:	Cy5
Immunogen Species:	Mitocryptide-1 antibody was raised against Human
Antigen Type:	Synthetic peptide
Immunogen:	Mitocryptide-1 antibody was raised against synthetic peptide corresponding to human Mitocryptide-1.
Specificity:	Human Mitocryptide-1.
Reactivity:	Human
Purification:	Purified
Presentation:	Liquid.
Recommended Storage:	Short term 4°C, long term aliquot and store at -20°C, avoid freeze thaw cycles. Sensitive to light.
Uses:	ELISA (Optimal dilution to be determined by the researcher)
Size:	100 µl
Requested From:	Japan

Laboratory Reagent For In Vitro Research Use Only

Not for resale without prior written consent from LifeSpan BioSciences, Inc.

Created on 9/27/2014

© 2014 LifeSpan BioSciences