

**Mitocryptide-1 Rabbit anti-Pig Polyclonal (Cy5) Antibody - LS-C203534 - LSBio**

<b>CatalogID:</b>	LS-C203534
<b>Target:</b>	Microcryptide-1
<b>Host</b>	Microcryptide-1 antibody was produced in Rabbit
<b>Clonality:</b>	Polyclonal
<b>Isotype:</b>	IgG
<b>Conjugations:</b>	Cy5
<b>Immunogen Species:</b>	Mitocryptide-1 antibody was raised against Pig
<b>Antigen Type:</b>	Synthetic peptide
<b>Immunogen:</b>	Mitocryptide-1 antibody was raised against synthetic peptide corresponding to porcine Mitocryptide-1.
<b>Specificity:</b>	Porcine Mitocryptide-1.
<b>Reactivity:</b>	Pig
<b>Purification:</b>	Purified
<b>Presentation:</b>	Liquid.
<b>Recommended Storage:</b>	Short term 4°C, long term aliquot and store at -20°C, avoid freeze thaw cycles. Sensitive to light.
<b>Uses:</b>	ELISA (Optimal dilution to be determined by the researcher)
<b>Size:</b>	100 µl
<b>Requested From:</b>	Japan

Laboratory Reagent For In Vitro Research Use Only

Not for resale without prior written consent from LifeSpan BioSciences, Inc.

Created on 9/27/2014

© 2014 LifeSpan BioSciences