

Hemagglutinin / HA Tag Mouse Monoclonal (FITC) Antibody - LS-C203621 - LSBio

CatalogID:	LS-C203621
Target:	HA Tag / Hemagglutinin Tag
Host	Mouse
Clonality:	Monoclonal
Isotype:	IgG1,k
Conjugations:	Fluorescein (FITC)
Antigen Type:	Synthetic peptide
Immunogen:	A synthetic peptide YPYDVPDYA coupled to KLH
Specificity:	C-terminal, N-terminal, and internal HA tagged fusion proteins.
Purification:	Protein A purified
Reconstitution:	Distilled water. Possible additional volumes for resuspension: 100 µl
Presentation:	Lyophilized from PBS, pH 7.4, 1% BSA, 0.02% sodium azide.
Recommended Storage:	Store lyophilized at -20°C. The reconstituted product can be stored for short term at 4 °C or long term at -20 °C or below. Avoid freeze/thaw cycles.
Uses:	Flow Cytometry (Optimal dilution to be determined by the researcher)
Size:	100 µg
Requested From:	Japan

Laboratory Reagent For In Vitro Research Use Only

Not for resale without prior written consent from LifeSpan BioSciences, Inc.

Created on 9/27/2014

© 2014 LifeSpan BioSciences