

CD20 Mouse anti-Human Monoclonal (FITC) (2H7) Antibody - LS-C204416 - LSBio

CatalogID:	LS-C204416
Target:	membrane-spanning 4-domains, subfamily A, member 1 (MS4A1)
Synonyms:	MS4A1 Antibody, B-lymphocyte antigen CD20 Antibody, B1 Antibody, Bp35 Antibody, CD20 Antibody, CD20 antigen Antibody, CD20 receptor Antibody, LEU-16 Antibody, CVID5 Antibody
Host	MS4A1 antibody was produced in Mouse
Clonality:	Monoclonal
Isotype:	IgG2b
Clone Name:	2H7
Conjugations:	Fluorescein (FITC)
Immunogen Species:	CD20 antibody was raised against Human
Antigen Type:	Cells
Immunogen:	CD20 antibody was raised against human tonsillar B cells.
Specificity:	Monoclonal antibody 2H7 recognizes CD20 (B1, Bp35), a 33-37 kD non-glycosylated membrane receptor with four transmembrane domains, expressed on pre-B lymphocytes, resting and activated B cells (not plasma cells), follicular dendritic cells, and at low levels on peripheral blood T lymphocytes.
Reactivity:	Human, Primate
Purification:	Protein A purified
Presentation:	PBS, 0.2% BSA, 0.02% sodium azide.
Recommended Storage:	Store in the dark at 2-8 C. Do not freeze. Avoid prolonged exposure to light.
Usage Summary:	The reagent is designed for Flow Cytometry analysis of human blood cells using 20 ul reagent / 100 ul of whole blood or 1 million cells in a suspension. The content of a vial (2 ml) is sufficient for 100 tests.
Uses:	Flow Cytometry (Optimal dilution to be determined by the researcher)
Size:	100 tst
Requested From:	Japan

Laboratory Reagent For In Vitro Research Use Only

Not for resale without prior written consent from LifeSpan BioSciences, Inc.

Created on 9/27/2014

© 2014 LifeSpan BioSciences