

CR1 / CD35 Mouse anti-Human Monoclonal (FITC) (E11) Antibody - LS-C204439 - LSBio	
CatalogID:	LS-C204439
Target:	complement component (3b/4b) receptor 1 (Knops blood group) (CR1)
Synonyms:	CR1 Antibody, C3b/C4b receptor Antibody, C3-binding protein Antibody, CD35 Antibody, C3BR Antibody, Knops blood group antigen Antibody, C4BR Antibody, CD35 antigen Antibody, Complement receptor type 1 Antibody, KN Antibody
Family / Subfamily:	Protease / Serine S1
Host	CR1 antibody was produced in Mouse
Clonality:	Monoclonal
Isotype:	IgG1
Clone Name:	E11
Conjugations:	Fluorescein (FITC)
Immunogen Species:	CR1 / CD35 antibody was raised against Human
Antigen Type:	Cells
Immunogen:	CR1 / CD35 antibody was raised against acute monocytic leukemia cells and normal blood monocytes.
Specificity:	Monoclonal antibody E11 recognizes CD35 (CR1), a type I glycoprotein expressed on granulocytes, monocytes, B cells, follicular dendritic cells, erythrocytes, NK and T cell subsets, as well as e.g. on glomerular podocytes.
Reactivity:	Human, Primate
Purification:	Protein A purified
Presentation:	PBS, 0.2% BSA, 0.02% sodium azide.
Recommended Storage:	Store in the dark at 2-8 C. Do not freeze. Avoid prolonged exposure to light.
Usage Summary:	The reagent is designed for Flow Cytometry analysis of human blood cells using 20 ul reagent / 100 ul of whole blood or 1 million cells in a suspension. The content of a vial (2 ml) is sufficient for 100 tests.
Uses:	Flow Cytometry (Optimal dilution to be determined by the researcher)
Size:	100 tst
Requested From:	Japan
Laboratory Reagent For In Vitro Research Use Only	
Not for resale without prior written consent from LifeSpan BioSciences, Inc.	
Created on 9/27/2014	
© 2014 LifeSpan BioSciences	