

SELP / P-Selectin / CD62P Mouse anti-Human Monoclonal (PE) (AK4) Antibody - LS-C204468 - LSBio	
CatalogID:	LS-C204468
Target:	selectin P (granule membrane protein 140kDa, antigen CD62) (SELP)
Synonyms:	SELP Antibody, CD62P Antibody, CD62 Antibody, GMP-140 Antibody, GMRP Antibody, Granule membrane protein 140 Antibody, LECAM3 Antibody, P-selectin Antibody, PSEL Antibody, CD62P antigen Antibody, GMP140 Antibody, Granulocyte membrane protein Antibody, GRMP Antibody, PADGEM Antibody
Host	SELP antibody was produced in Mouse
Clonality:	Monoclonal
Isotype:	IgG1
Clone Name:	AK4
Conjugations:	Phycoerythrin (PE)
Immunogen Species:	SELP / P-Selectin / CD62P antibody was raised against Human
Antigen Type:	Cells
Immunogen:	SELP / P-Selectin / CD62P antibody was raised against human platelets.
Specificity:	Monoclonal antibody HI62P recognizes CD62P (P-selectin), a 140 kD single chain type I transmembrane glycoprotein present in secretory alpha-granules in platelets, in Weibel-Palade bodies in endothelial cells and in megakaryocytes; it is relocated to the plasma membrane upon activation.
Reactivity:	Human, Primate
Purification:	Protein A purified
Presentation:	PBS, 0.2% BSA, 0.02% sodium azide.
Recommended Storage:	Store in the dark at 2-8 C. Do not freeze. Avoid prolonged exposure to light.
Usage Summary:	The reagent is designed for Flow Cytometry analysis of human blood cells using 20 ul reagent / 100 ul of whole blood or 1 million cells in a suspension. The content of a vial (2 ml) is sufficient for 100 tests.
Uses:	Flow Cytometry (Optimal dilution to be determined by the researcher)
Size:	100 tst
Requested From:	Japan
Laboratory Reagent For In Vitro Research Use Only	
Not for resale without prior written consent from LifeSpan BioSciences, Inc.	
Created on 9/27/2014	
© 2014 LifeSpan BioSciences	