

BCL2 / Bcl-2 Mouse anti-Human Monoclonal (aa41-54) (PE) (Bcl2/100) Antibody - LS-C204472 - LSBio	
CatalogID:	LS-C204472
Target:	B-cell CLL/lymphoma 2 (BCL2)
Synonyms:	BCL2 Antibody, Apoptosis regulator Bcl-2 Antibody, B-cell CLL/lymphoma 2 Antibody, B-cell lymphoma protein 2 Antibody, Bcl-2 Antibody, PPP1R50 Antibody
Family / Subfamily:	Apoptosis
Host	BCL2 antibody was produced in Mouse
Clonality:	Monoclonal
Isotype:	IgG1
Clone Name:	Bcl2/100
Conjugations:	Phycoerythrin (PE)
Immunogen Species:	BCL2 / Bcl-2 antibody was raised against Human
Antigen Type:	Synthetic peptide
Immunogen:	BCL2 / Bcl-2 antibody was raised against synthetic peptide corresponding to the amino acids 41-54 of human Bcl2.
Specificity:	Monoclonal antibody Bcl2/100 recognizes Bcl2, a 26 kD protooncogene with anti-apoptotic effect, expressed in outer mitochondrial membrane, endoplasmic reticulum and nuclear envelope.
Epitope:	aa41-54
Reactivity:	Human
Purification:	Protein A purified
Presentation:	PBS, 0.2% BSA, 0.02% sodium azide.
Recommended Storage:	Store in the dark at 2-8 C. Do not freeze. Avoid prolonged exposure to light.
Usage Summary:	The reagent is designed for Flow Cytometry analysis of human blood cells using 20 ul reagent / 100 ul of whole blood or 1 million cells in a suspension. The content of a vial (2 ml) is sufficient for 100 tests.
Uses:	Flow Cytometry (Optimal dilution to be determined by the researcher)
Size:	100 tst
Requested From:	Japan
Laboratory Reagent For In Vitro Research Use Only	
Not for resale without prior written consent from LifeSpan BioSciences, Inc.	
Created on 9/27/2014	
© 2014 LifeSpan BioSciences	