

HLA-DR Mouse anti-Human Monoclonal (PE) (L243) Antibody - LS-C204477 - LSBio

CatalogID:	LS-C204477
Target:	MHC Class II HLA-DR
Host	MHC Class II HLA-DR antibody was produced in Mouse
Clonality:	Monoclonal
Isotype:	IgG2a
Clone Name:	L243
Conjugations:	Phycoerythrin (PE)
Immunogen Species:	HLA-DR antibody was raised against Human
Antigen Type:	Cells
Immunogen:	HLA-DR antibody was raised against human B lymphocytes.
Specificity:	Monoclonal antibody L243 recognizes specifically HLA-DR molecules, both peptide-loaded and empty.
Reactivity:	Human, Dog, Primate
Purification:	Protein A purified
Presentation:	PBS, 0.2% BSA, 0.02% sodium azide.
Recommended Storage:	Store in the dark at 2-8 C. Do not freeze. Avoid prolonged exposure to light.
Usage Summary:	The reagent is designed for Flow Cytometry analysis of human blood cells using 10 ul reagent / 100 ul of whole blood or 1 million cells in a suspension. The content of a vial (1 ml) is sufficient for 100 tests.
Uses:	Flow Cytometry (Optimal dilution to be determined by the researcher)
Size:	100 tst
Requested From:	Japan

Laboratory Reagent For In Vitro Research Use Only

Not for resale without prior written consent from LifeSpan BioSciences, Inc.

Created on 9/27/2014

© 2014 LifeSpan BioSciences