

CD5 Mouse anti-Human Monoclonal (PE) (L17F12) Antibody - LS-C204483 - LSBio

<b>CatalogID:</b>	LS-C204483
<b>Target:</b>	CD5 molecule
<b>Synonyms:</b>	CD5 Antibody, CD5 antigen Antibody, CD5 antigen (p56-62) Antibody, CD5 molecule Antibody, Lymphocyte antigen T1/Leu-1 Antibody, LEU1 Antibody, T1 Antibody
<b>Host</b>	CD5 antibody was produced in Mouse
<b>Clonality:</b>	Monoclonal
<b>Isotype:</b>	IgG2a
<b>Clone Name:</b>	L17F12
<b>Conjugations:</b>	Phycoerythrin (PE)
<b>Immunogen Species:</b>	CD5 antibody was raised against Human
<b>Antigen Type:</b>	Cells
<b>Immunogen:</b>	CD5 antibody was raised against human acute lymphoblastic leukemia (ALL) T cells.
<b>Specificity:</b>	Monoclonal antibody L17F12 reacts with CD5, a 67kD single-chain transmembrane glycoprotein expressed on mature T lymphocytes, most of thymocytes and B lymphocytes subset (B-1a lymphocytes).
<b>Reactivity:</b>	Human
<b>Purification:</b>	Protein A purified
<b>Presentation:</b>	PBS, 0.2% BSA, 0.02% sodium azide.
<b>Recommended Storage:</b>	Store in the dark at 2-8 C. Do not freeze. Avoid prolonged exposure to light.
<b>Usage Summary:</b>	The reagent is designed for Flow Cytometry analysis of human blood cells using 10 ul reagent / 100 ul of whole blood or 1 million cells in a suspension. The content of a vial (1 ml) is sufficient for 100 tests.
<b>Uses:</b>	Flow Cytometry (Optimal dilution to be determined by the researcher)
<b>Size:</b>	100 tst
<b>Requested From:</b>	Japan

Laboratory Reagent For In Vitro Research Use Only

Not for resale without prior written consent from LifeSpan BioSciences, Inc.

Created on 9/27/2014

© 2014 LifeSpan BioSciences