

CR2 / CD21 Rabbit anti-Human Polyclonal (C-Terminus) (FITC) Antibody - LS-C206425 - LSBio	
CatalogID:	LS-C206425
Target:	complement component (3d/Epstein Barr virus) receptor 2 (CR2)
Synonyms:	CR2 Antibody, C3DR Antibody, CD21 antigen Antibody, CVID7 Antibody, Complement receptors type 2 Antibody, CD21 Antibody, Complement C3d receptor Antibody, EBV receptor Antibody, SLEB9 Antibody, Complement receptor type 2 Antibody, CR Antibody, Epstein-Barr virus receptor Antibody
Host	CR2 antibody was produced in Rabbit
Clonality:	Polyclonal
Isotype:	IgG
Conjugations:	Fluorescein (FITC)
Immunogen Species:	CR2 / CD21 antibody was raised against Human
Antigen Type:	Synthetic peptide
Immunogen:	CR2 / CD21 antibody was raised against kLH-conjugated synthetic peptide from the C-terminus of human CR2.
Specificity:	Human CR2 / CD21
Epitope:	C-Terminus
Reactivity:	Human
Purification:	Immunoaffinity purified
Presentation:	PBS, pH 7.2, 0.09% sodium azide.
Recommended Storage:	May be stored at 4°C for short-term only. Aliquot to avoid repeated freezing and thawing. Store at -20°C. Aliquots are stable for at least 6 months.
Usage Summary:	The applications listed above are for the unconjugated form of this antibody. The conjugated antibody has not been tested.
Uses:	IHC, Western blot, Flow Cytometry, ELISA (Optimal dilution to be determined by the researcher)
Size:	200 µl
Requested From:	Japan
Laboratory Reagent For In Vitro Research Use Only	
Not for resale without prior written consent from LifeSpan BioSciences, Inc.	
Created on 9/27/2014	
© 2014 LifeSpan BioSciences	