

CR2 / CD21 Rabbit anti-Human Polyclonal (C-Terminus) (PE) Antibody - LS-C206430 - LSBio

<b>CatalogID:</b>	LS-C206430
<b>Target:</b>	complement component (3d/Epstein Barr virus) receptor 2 (CR2)
<b>Synonyms:</b>	CR2 Antibody, C3DR Antibody, CD21 antigen Antibody, CVID7 Antibody, Complement receptors type 2 Antibody, CD21 Antibody, Complement C3d receptor Antibody, EBV receptor Antibody, SLEB9 Antibody, Complement receptor type 2 Antibody, CR Antibody, Epstein-Barr virus receptor Antibody
<b>Host</b>	CR2 antibody was produced in Rabbit
<b>Clonality:</b>	Polyclonal
<b>Isotype:</b>	IgG
<b>Conjugations:</b>	Phycoerythrin (PE)
<b>Immunogen Species:</b>	CR2 / CD21 antibody was raised against Human
<b>Antigen Type:</b>	Synthetic peptide
<b>Immunogen:</b>	CR2 / CD21 antibody was raised against kLH-conjugated synthetic peptide from the C-terminus of human CR2.
<b>Specificity:</b>	Human CR2 / CD21
<b>Epitope:</b>	C-Terminus
<b>Reactivity:</b>	Human
<b>Purification:</b>	Immunoaffinity purified
<b>Presentation:</b>	PBS, pH 7.2, 0.09% sodium azide.
<b>Recommended Storage:</b>	May be stored at 4°C for short-term only. Aliquot to avoid repeated freezing and thawing. Store at -20°C. Aliquots are stable for at least 6 months.
<b>Usage Summary:</b>	The applications listed above are for the unconjugated form of this antibody. The conjugated antibody has not been tested.
<b>Uses:</b>	IHC, Western blot, Flow Cytometry, ELISA (Optimal dilution to be determined by the researcher)
<b>Size:</b>	200 µl
<b>Requested From:</b>	Japan

Laboratory Reagent For In Vitro Research Use Only

Not for resale without prior written consent from LifeSpan BioSciences, Inc.

Created on 9/27/2014

© 2014 LifeSpan BioSciences