

CKM / Creatine Kinase MM Rabbit anti-Human Polyclonal (C-Terminus) (FITC) Antibody - LS-C206549 - LSBio

CatalogID:	LS-C206549
Target:	creatine kinase, muscle (CKM)
Synonyms:	CKM Antibody, CK-MB Antibody, CKMM Antibody, Creatine kinase M-type Antibody, Creatine kinase, muscle Antibody, M-CK Antibody, Creatine kinase M chain Antibody
Host	CKM antibody was produced in Rabbit
Clonality:	Polyclonal
Isotype:	IgG
Conjugations:	Fluorescein (FITC)
Immunogen Species:	CKM / Creatine Kinase MM antibody was raised against Human
Antigen Type:	Synthetic peptide
Immunogen:	CKM / Creatine Kinase MM antibody was raised against kLH-conjugated synthetic peptide from the C-terminus of human CKM.
Specificity:	Human CKM / Creatine Kinase MM
Epitope:	C-Terminus
Reactivity:	Human
Purification:	Protein G purified
Presentation:	PBS, pH 7.2, 0.09% sodium azide.
Recommended Storage:	May be stored at 4°C for short-term only. Aliquot to avoid repeated freezing and thawing. Store at -20°C. Aliquots are stable for at least 6 months.
Usage Summary:	The applications listed above are for the unconjugated form of this antibody. The conjugated antibody has not been tested.
Uses:	Western blot, ELISA (Optimal dilution to be determined by the researcher)
Size:	200 µl
Requested From:	Japan

Laboratory Reagent For In Vitro Research Use Only

Not for resale without prior written consent from LifeSpan BioSciences, Inc.

Created on 9/27/2014

© 2014 LifeSpan BioSciences