

FAU Rabbit anti-Human Polyclonal (Internal) (FITC) Antibody - LS-C206616 - LSBio

CatalogID:	LS-C206616
Target:	Finkel-Biskis-Reilly murine sarcoma virus (FBR-MuSV) ubiquitously expressed (FAU)
Synonyms:	FAU Antibody, 40S ribosomal protein S30 Antibody, Asr1 Antibody, Fub1 Antibody, MNSFbeta Antibody, S30 Antibody, FAU1 Antibody, Fubi Antibody, Ribosomal protein S30 Antibody, RPS30 Antibody
Host	FAU antibody was produced in Rabbit
Clonality:	Polyclonal
Isotype:	IgG
Conjugations:	Fluorescein (FITC)
Immunogen Species:	FAU antibody was raised against Human
Antigen Type:	Synthetic peptide
Immunogen:	FAU antibody was raised against kLH-conjugated synthetic peptide from internal region of human FUBI.
Specificity:	Human FAU
Epitope:	Internal
Reactivity:	Human
Purification:	Protein G purified
Presentation:	PBS, pH 7.2, 0.09% sodium azide.
Recommended Storage:	May be stored at 4°C for short-term only. Aliquot to avoid repeated freezing and thawing. Store at -20°C. Aliquots are stable for at least 6 months.
Usage Summary:	The applications listed above are for the unconjugated form of this antibody. The conjugated antibody has not been tested.
Uses:	IHC, Western blot, ELISA (Optimal dilution to be determined by the researcher)
Size:	200 µl
Requested From:	Japan

Laboratory Reagent For In Vitro Research Use Only

Not for resale without prior written consent from LifeSpan BioSciences, Inc.

Created on 9/27/2014

© 2014 LifeSpan BioSciences