

FOS / c-FOS Rabbit anti-Human Polyclonal (pThr232) (PE) Antibody - LS-C206635 - LSBio	
CatalogID:	LS-C206635
Target:	FBJ murine osteosarcoma viral oncogene homolog (FOS)
Synonyms:	FOS Antibody, Activator protein 1 Antibody, AP-1 Antibody, Cellular oncogene c-fos Antibody, Cellular oncogene fos Antibody, Proto-oncogene c-Fos Antibody, v-fos Antibody, C-FOS Antibody, G0s7 Antibody, p55 Antibody
Host	FOS antibody was produced in Rabbit
Clonality:	Polyclonal
Isotype:	IgG
Conjugations:	Phycoerythrin (PE)
Immunogen Species:	FOS / c-FOS antibody was raised against Human
Antigen Type:	Synthetic phospho-peptide
Immunogen:	FOS / c-FOS antibody was raised against kLH-conjugated synthetic phosphopeptide around Thr232 of human FOS.
Specificity:	Human FOS / c-FOS
Epitope:	pThr232
Reactivity:	Human
Purification:	Protein A purified
Presentation:	PBS, pH 7.2, 0.09% sodium azide.
Recommended Storage:	May be stored at 4°C for short-term only. Aliquot to avoid repeated freezing and thawing. Store at -20°C. Aliquots are stable for at least 6 months.
Usage Summary:	The applications listed above are for the unconjugated form of this antibody. The conjugated antibody has not been tested.
Uses:	ELISA, Dot Blot (Optimal dilution to be determined by the researcher)
Size:	200 µl
Requested From:	Japan
Laboratory Reagent For In Vitro Research Use Only	
Not for resale without prior written consent from LifeSpan BioSciences, Inc.	
Created on 9/27/2014	
© 2014 LifeSpan BioSciences	