

CDKN1A / WAF1 / p21 Rabbit anti-Human Polyclonal (pThr145) (FITC) Antibody - LS-C206780 - LSBio	
<b>CatalogID:</b>	LS-C206780
<b>Target:</b>	cyclin-dependent kinase inhibitor 1A (p21, Cip1) (CDKN1A)
<b>Synonyms:</b>	CDKN1A Antibody, CDKN1 Antibody, CIP1 Antibody, CDK-interacting protein 1 Antibody, DNA synthesis inhibitor Antibody, p21CIP1 Antibody, MDA6 Antibody, PIC1 Antibody, SDI1 Antibody, WAF1 Antibody, CAP20 Antibody, CDK-interaction protein 1 Antibody, MDA-6 Antibody, p21 Antibody
<b>Host</b>	CDKN1A antibody was produced in Rabbit
<b>Clonality:</b>	Polyclonal
<b>Isotype:</b>	IgG
<b>Conjugations:</b>	Fluorescein (FITC)
<b>Immunogen Species:</b>	CDKN1A / WAF1 / p21 antibody was raised against Human
<b>Antigen Type:</b>	Synthetic phospho-peptide
<b>Immunogen:</b>	CDKN1A / WAF1 / p21 antibody was raised against kLH-conjugated synthetic phosphopeptide around Thr145 of human P21CIP1.
<b>Specificity:</b>	Human CDKN1A / WAF1 / p21
<b>Epitope:</b>	pThr145
<b>Reactivity:</b>	Human
<b>Purification:</b>	Immunoaffinity purified
<b>Presentation:</b>	PBS, pH 7.2, 0.09% sodium azide.
<b>Recommended Storage:</b>	May be stored at 4°C for short-term only. Aliquot to avoid repeated freezing and thawing. Store at -20°C. Aliquots are stable for at least 6 months.
<b>Usage Summary:</b>	The applications listed above are for the unconjugated form of this antibody. The conjugated antibody has not been tested.
<b>Uses:</b>	IHC, Western blot, ELISA, Dot Blot (Optimal dilution to be determined by the researcher)
<b>Size:</b>	200 µl
<b>Requested From:</b>	Japan
Laboratory Reagent For In Vitro Research Use Only	
Not for resale without prior written consent from LifeSpan BioSciences, Inc.	
Created on 9/27/2014	
© 2014 LifeSpan BioSciences	