

CDKN1B / p27 Kip1 Rabbit anti-Human Polyclonal (pSer10) (FITC) Antibody - LS-C206794 - LSBio	
CatalogID:	LS-C206794
Target:	cyclin-dependent kinase inhibitor 1B (p27, Kip1) (CDKN1B)
Synonyms:	CDKN1B Antibody, KIP1 Antibody, MEN1B Antibody, P27 Kip1 Antibody, p27KIP1 Antibody, MEN4 Antibody, CDKN4 Antibody
Host	CDKN1B antibody was produced in Rabbit
Clonality:	Polyclonal
Isotype:	IgG
Conjugations:	Fluorescein (FITC)
Immunogen Species:	CDKN1B / p27 Kip1 antibody was raised against Human
Antigen Type:	Synthetic phospho-peptide
Immunogen:	CDKN1B / p27 Kip1 antibody was raised against kLH-conjugated synthetic phosphopeptide around Ser10 of human p27Kip1.
Specificity:	Human CDKN1B / p27 Kip1
Epitope:	pSer10
Reactivity:	Human
Purification:	Immunoaffinity purified
Presentation:	PBS, pH 7.2, 0.09% sodium azide.
Recommended Storage:	May be stored at 4°C for short-term only. Aliquot to avoid repeated freezing and thawing. Store at -20°C. Aliquots are stable for at least 6 months.
Usage Summary:	The applications listed above are for the unconjugated form of this antibody. The conjugated antibody has not been tested.
Uses:	IHC, Western blot, ELISA (Optimal dilution to be determined by the researcher)
Size:	200 µl
Requested From:	Japan
Laboratory Reagent For In Vitro Research Use Only	
Not for resale without prior written consent from LifeSpan BioSciences, Inc.	
Created on 9/27/2014	
© 2014 LifeSpan BioSciences	