

CDKN1B / p27 Kip1 Rabbit anti-Human Polyclonal (pThr198) (FITC) Antibody - LS-C206813 - LSBio

CatalogID:	LS-C206813
Target:	cyclin-dependent kinase inhibitor 1B (p27, Kip1) (CDKN1B)
Synonyms:	CDKN1B Antibody, KIP1 Antibody, MEN1B Antibody, P27 Kip1 Antibody, p27KIP1 Antibody, MEN4 Antibody, CDKN4 Antibody
Host	CDKN1B antibody was produced in Rabbit
Clonality:	Polyclonal
Isotype:	IgG
Conjugations:	Fluorescein (FITC)
Immunogen Species:	CDKN1B / p27 Kip1 antibody was raised against Human
Antigen Type:	Synthetic phospho-peptide
Immunogen:	CDKN1B / p27 Kip1 antibody was raised against kLH-conjugated synthetic phosphopeptide around Thr198 of human 27Kip1.
Specificity:	Human CDKN1B / p27 Kip1
Epitope:	pThr198
Reactivity:	Human
Purification:	Protein A purified
Presentation:	PBS, pH 7.2, 0.09% sodium azide.
Recommended Storage:	May be stored at 4°C for short-term only. Aliquot to avoid repeated freezing and thawing. Store at -20°C. Aliquots are stable for at least 6 months.
Usage Summary:	The applications listed above are for the unconjugated form of this antibody. The conjugated antibody has not been tested.
Uses:	ELISA, Dot Blot (Optimal dilution to be determined by the researcher)
Size:	200 µl
Requested From:	Japan

Laboratory Reagent For In Vitro Research Use Only

Not for resale without prior written consent from LifeSpan BioSciences, Inc.

Created on 9/27/2014

© 2014 LifeSpan BioSciences