

CDKN2A / p16INK4a Rabbit anti-Human Polyclonal (aa112-141) (FITC) Antibody - LS-C207245 - LSBio

CatalogID:	LS-C207245
Target:	cyclin-dependent kinase inhibitor 2A (CDKN2A)
Synonyms:	CDKN2A Antibody, ARF Antibody, Cdk4 inhibitor Antibody, CDK4 inhibitor p16-INK4 Antibody, INK4 Antibody, INK4A Antibody, Multiple tumor suppressor 1 Antibody, MTS-1 Antibody, MTS1 Antibody, p16-INK4 Antibody, p16-INK4A Antibody, p16INK4 Antibody, p19 Antibody, p14 Antibody, p16 Antibody, p16 INK4A Antibody, p16INK4A Antibody, p14ARF Antibody, TP16 Antibody, CDK4I Antibody, CDKN2 Antibody, CMM2 Antibody, MLM Antibody, p19ARF Antibody
Host	CDKN2A antibody was produced in Rabbit
Clonality:	Polyclonal
Isotype:	IgG
Conjugations:	Fluorescein (FITC)
Immunogen Species:	CDKN2A / p16INK4a antibody was raised against Human
Antigen Type:	Synthetic peptide
Immunogen:	CDKN2A / p16INK4a antibody was raised against kLH-conjugated synthetic peptide (aa112-141) from the C-terminus of human MLM.
Specificity:	Human CDKN2A / p14ARF
Epitope:	aa112-141
Reactivity:	Human
Purification:	Protein A purified
Presentation:	PBS, pH 7.2, 0.09% sodium azide.
Recommended Storage:	May be stored at 4°C for short-term only. Aliquot to avoid repeated freezing and thawing. Store at -20°C. Aliquots are stable for 6 months.
Usage Summary:	The applications listed above are for the unconjugated form of this antibody. The conjugated antibody has not been tested.
Uses:	IHC, Western blot, Flow Cytometry, ELISA (Optimal dilution to be determined by the researcher)
Size:	200 µl
Requested From:	Japan

Laboratory Reagent For In Vitro Research Use Only

Not for resale without prior written consent from LifeSpan BioSciences, Inc.

Created on 9/27/2014

© 2014 LifeSpan BioSciences