

**CDKN1A / WAF1 / p21 Rabbit anti-Mouse Polyclonal (pSer78) (FITC) Antibody - LS-C207394 - LSBio**

<b>CatalogID:</b>	LS-C207394
<b>Target:</b>	cyclin-dependent kinase inhibitor 1A (p21, Cip1) (CDKN1A)
<b>Synonyms:</b>	CDKN1A Antibody, CDKN1 Antibody, CIP1 Antibody, CDK-interacting protein 1 Antibody, DNA synthesis inhibitor Antibody, p21CIP1 Antibody, MDA6 Antibody, PIC1 Antibody, SDI1 Antibody, WAF1 Antibody, CAP20 Antibody, CDK-interaction protein 1 Antibody, MDA-6 Antibody, p21 Antibody
<b>Host</b>	CDKN1A antibody was produced in Rabbit
<b>Clonality:</b>	Polyclonal
<b>Isotype:</b>	IgG
<b>Conjugations:</b>	Fluorescein (FITC)
<b>Immunogen Species:</b>	CDKN1A / WAF1 / p21 antibody was raised against Mouse
<b>Antigen Type:</b>	Synthetic phospho-peptide
<b>Immunogen:</b>	CDKN1A / WAF1 / p21 antibody was raised against kLH-conjugated synthetic phosphopeptide around Ser78 of mouse p21Cip1.
<b>Specificity:</b>	Mouse CDKN1A / WAF1 / p21
<b>Epitope:</b>	pSer78
<b>Reactivity:</b>	Mouse
<b>Purification:</b>	Protein A purified
<b>Presentation:</b>	PBS, pH 7.2, 0.09% sodium azide.
<b>Recommended Storage:</b>	May be stored at 4°C for short-term only. Aliquot to avoid repeated freezing and thawing. Store at -20°C. Aliquots are stable for 6 months.
<b>Usage Summary:</b>	The applications listed above are for the unconjugated form of this antibody. The conjugated antibody has not been tested.
<b>Uses:</b>	ELISA, Dot Blot (Optimal dilution to be determined by the researcher)
<b>Size:</b>	200 µl
<b>Requested From:</b>	Japan

Laboratory Reagent For In Vitro Research Use Only

Not for resale without prior written consent from LifeSpan BioSciences, Inc.

Created on 9/27/2014

© 2014 LifeSpan BioSciences