

SERPINA6 / CBG Rabbit anti-Human Polyclonal (aa200-229) (PE) Antibody - LS-C207638 - LSBio	
CatalogID:	LS-C207638
Target:	serpin peptidase inhibitor, clade A (alpha-1 antiproteinase, antitrypsin), member 6 (SERPINA6)
Synonyms:	SERPINA6 Antibody, Transcortin Antibody, CBG Antibody, Serpin A6 Antibody
Host	SERPINA6 antibody was produced in Rabbit
Clonality:	Polyclonal
Isotype:	IgG
Conjugations:	Phycoerythrin (PE)
Immunogen Species:	SERPINA6 / CBG antibody was raised against Human
Antigen Type:	Synthetic peptide
Immunogen:	SERPINA6 / CBG antibody was raised against KLH-conjugated synthetic peptide (aa200-229) from human SERPINA6.
Specificity:	Human SERPINA6 / CBG.
Epitope:	aa200-229
Reactivity:	Human, Mouse
Purification:	Protein A purified
Presentation:	PBS, pH 7.2, 0.09% sodium azide.
Recommended Storage:	May be stored at 4°C for short-term only. Aliquot to avoid repeated freezing and thawing. Store at -20°C. Aliquots are stable for 6 months.
Usage Summary:	The applications listed above are for the unconjugated form of this antibody. The conjugated antibody has not been tested.
Uses:	Immunofluorescence, Western blot, ELISA (Optimal dilution to be determined by the researcher)
Size:	200 µl
Requested From:	Japan
Laboratory Reagent For In Vitro Research Use Only	
Not for resale without prior written consent from LifeSpan BioSciences, Inc.	
Created on 9/27/2014	
© 2014 LifeSpan BioSciences	