

VRL1 / TRPV2 Rabbit anti-Human Polyclonal (aa51-81) (PE) Antibody - LS-C207820 - LSBio	
CatalogID:	LS-C207820
Target:	transient receptor potential cation channel, subfamily V, member 2 (TRPV2)
Synonyms:	TRPV2 Antibody, GRC Antibody, OTRPC2 Antibody, Osm-9-like TRP channel 2 Antibody, VRL-1 Antibody, VRL1 Antibody, VRL Antibody
Family / Subfamily:	Ion Channel / Calcium channel - TRP
Host	TRPV2 antibody was produced in Rabbit
Clonality:	Polyclonal
Isotype:	IgG
Conjugations:	Phycoerythrin (PE)
Immunogen Species:	VRL1 / TRPV2 antibody was raised against Human
Antigen Type:	Synthetic peptide
Immunogen:	VRL1 / TRPV2 antibody was raised against kLH-conjugated synthetic peptide (aa51-81) from the N-terminus of human TRPV2.
Specificity:	Human VRL1 / TRPV2
Epitope:	aa51-81
Reactivity:	Human
Purification:	Protein A purified
Presentation:	PBS, pH 7.2, 0.09% sodium azide.
Recommended Storage:	May be stored at 4°C for short-term only. Aliquot to avoid repeated freezing and thawing. Store at -20°C. Aliquots are stable for 6 months.
Usage Summary:	The applications listed above are for the unconjugated form of this antibody. The conjugated antibody has not been tested.
Uses:	Western blot, ELISA (Optimal dilution to be determined by the researcher)
Size:	200 µl
Requested From:	Japan
Laboratory Reagent For In Vitro Research Use Only	
Not for resale without prior written consent from LifeSpan BioSciences, Inc.	
Created on 9/27/2014	
© 2014 LifeSpan BioSciences	