

VEGFA / VEGF Rabbit anti-Human Polyclonal (aa73-103) (PE) Antibody - LS-C207852 - LSBio

CatalogID:	LS-C207852
Target:	vascular endothelial growth factor A (VEGFA)
Synonyms:	VEGFA Antibody, VPF Antibody, Vascular permeability factor Antibody, VEGF Antibody, VEGF-A Antibody, MVCD1 Antibody
Family / Subfamily:	PDGF / not assigned-PDGF
Host	VEGFA antibody was produced in Rabbit
Clonality:	Polyclonal
Isotype:	IgG
Conjugations:	Phycoerythrin (PE)
Immunogen Species:	VEGFA / VEGF antibody was raised against Human
Antigen Type:	Synthetic peptide
Immunogen:	VEGFA / VEGF antibody was raised against kLH-conjugated synthetic peptide (aa73-103) from the N-terminus of human VEGF1.
Specificity:	Human VEGFA / VEGF
Epitope:	aa73-103
Reactivity:	Human
Purification:	Protein A purified
Presentation:	PBS, pH 7.2, 0.09% sodium azide.
Recommended Storage:	May be stored at 4°C for short-term only. Aliquot to avoid repeated freezing and thawing. Store at -20°C. Aliquots are stable for 6 months.
Usage Summary:	The applications listed above are for the unconjugated form of this antibody. The conjugated antibody has not been tested.
Uses:	Western blot, ELISA (Optimal dilution to be determined by the researcher)
Size:	200 µl
Requested From:	Japan

Laboratory Reagent For In Vitro Research Use Only

Not for resale without prior written consent from LifeSpan BioSciences, Inc.

Created on 9/27/2014

© 2014 LifeSpan BioSciences