

VEGFA / VEGF Rabbit anti-Human Polyclonal (aa147-176) (FITC) Antibody - LS-C207860 - LSBio	
CatalogID:	LS-C207860
Target:	vascular endothelial growth factor A (VEGFA)
Synonyms:	VEGFA Antibody, VPF Antibody, Vascular permeability factor Antibody, VEGF Antibody, VEGF-A Antibody, MVCD1 Antibody
Family / Subfamily:	PDGF / not assigned-PDGF
Host	VEGFA antibody was produced in Rabbit
Clonality:	Polyclonal
Isotype:	IgG
Conjugations:	Fluorescein (FITC)
Immunogen Species:	VEGFA / VEGF antibody was raised against Human
Antigen Type:	Synthetic peptide
Immunogen:	VEGFA / VEGF antibody was raised against kLH-conjugated synthetic peptide (aa147-176) from human VEGFA.
Specificity:	Human VEGFA / VEGF
Epitope:	aa147-176
Reactivity:	Human
Purification:	Protein A purified
Presentation:	PBS, pH 7.2, 0.09% sodium azide.
Recommended Storage:	May be stored at 4°C for short-term only. Aliquot to avoid repeated freezing and thawing. Store at -20°C. Aliquots are stable for 6 months.
Usage Summary:	The applications listed above are for the unconjugated form of this antibody. The conjugated antibody has not been tested.
Uses:	Western blot, ELISA (Optimal dilution to be determined by the researcher)
Size:	200 µl
Requested From:	Japan
Laboratory Reagent For In Vitro Research Use Only	
Not for resale without prior written consent from LifeSpan BioSciences, Inc.	
Created on 9/27/2014	
© 2014 LifeSpan BioSciences	