

ATF4 Rabbit anti-Human Polyclonal (pSer245) (PE) Antibody - LS-C208076 - LSBio

CatalogID:	LS-C208076
Target:	activating transcription factor 4 (ATF4)
Synonyms:	ATF4 Antibody, CREB-2 Antibody, CREB2 Antibody, TXREB Antibody, DNA-binding protein TAXREB67 Antibody, TAXREB67 Antibody
Family / Subfamily:	Transcription factor / not assigned-Transcription factor
Host	ATF4 antibody was produced in Rabbit
Clonality:	Polyclonal
Isotype:	IgG
Conjugations:	Phycoerythrin (PE)
Immunogen Species:	ATF4 antibody was raised against Human
Antigen Type:	Synthetic phospho-peptide
Immunogen:	ATF4 antibody was raised against kLH-conjugated synthetic phosphopeptide around Ser245 of human ATF4.
Specificity:	Human ATF4
Epitope:	pSer245
Reactivity:	Human
Purification:	Immunoaffinity purified
Presentation:	PBS, pH 7.2, 0.09% sodium azide.
Recommended Storage:	May be stored at 4°C for short-term only. Aliquot to avoid repeated freezing and thawing. Store at -20°C. Aliquots are stable for 6 months.
Usage Summary:	The applications listed above are for the unconjugated form of this antibody. The conjugated antibody has not been tested.
Uses:	ELISA, Dot Blot (Optimal dilution to be determined by the researcher)
Size:	200 µl
Requested From:	Japan

Laboratory Reagent For In Vitro Research Use Only

Not for resale without prior written consent from LifeSpan BioSciences, Inc.

Created on 9/27/2014

© 2014 LifeSpan BioSciences