

BECN1 / Beclin-1 Rabbit anti-Human Polyclonal (aa106-136) (FITC) Antibody - LS-C208109 - LSBio

CatalogID:	LS-C208109
Target:	beclin 1, autophagy related (BECN1)
Synonyms:	BECN1 Antibody, ATG6 Antibody, Beclin 1 Antibody, Beclin 1, autophagy related Antibody, Beclin1 Antibody, GT197 Antibody, Protein GT197 Antibody, VPS30 Antibody, Beclin-1 Antibody
Host	BECN1 antibody was produced in Rabbit
Clonality:	Polyclonal
Isotype:	IgG
Conjugations:	Fluorescein (FITC)
Immunogen Species:	BECN1 / Beclin-1 antibody was raised against Human
Antigen Type:	Synthetic peptide
Immunogen:	BECN1 / Beclin-1 antibody was raised against kLH-conjugated synthetic peptide (aa106-136) from the human Beclin1-BH3.
Specificity:	Human BECN1 / Beclin-1
Epitope:	aa106-136
Reactivity:	Human, Mouse
Purification:	Protein A purified
Presentation:	PBS, pH 7.2, 0.09% sodium azide.
Recommended Storage:	May be stored at 4°C for short-term only. Aliquot to avoid repeated freezing and thawing. Store at -20°C. Aliquots are stable for 6 months.
Usage Summary:	The applications listed above are for the unconjugated form of this antibody. The conjugated antibody has not been tested.
Uses:	IHC, Western blot, Flow Cytometry, ELISA (Optimal dilution to be determined by the researcher)
Size:	200 µl
Requested From:	Japan

Laboratory Reagent For In Vitro Research Use Only

Not for resale without prior written consent from LifeSpan BioSciences, Inc.

Created on 9/27/2014

© 2014 LifeSpan BioSciences