

**CD55 Rabbit anti-Human Polyclonal (aa50-79) (FITC) Antibody - LS-C208288 - LSBio**

<b>CatalogID:</b>	LS-C208288
<b>Target:</b>	CD55 molecule, decay accelerating factor for complement (Cromer blood group)
<b>Synonyms:</b>	CD55 Antibody, CD55 antigen Antibody, Decay-accelerating factor Antibody, TC Antibody, CR Antibody, CROM Antibody, DAF Antibody
<b>Host</b>	CD55 antibody was produced in Rabbit
<b>Clonality:</b>	Polyclonal
<b>Isotype:</b>	IgG
<b>Conjugations:</b>	Fluorescein (FITC)
<b>Immunogen Species:</b>	CD55 antibody was raised against Human
<b>Antigen Type:</b>	Synthetic peptide
<b>Immunogen:</b>	CD55 antibody was raised against kLH-conjugated synthetic peptide (aa50-79) from the N-terminus of human CD55.
<b>Specificity:</b>	Human CD55
<b>Epitope:</b>	aa50-79
<b>Reactivity:</b>	Human
<b>Purification:</b>	Protein A purified
<b>Presentation:</b>	PBS, pH 7.2, 0.09% sodium azide.
<b>Recommended Storage:</b>	May be stored at 4°C for short-term only. Aliquot to avoid repeated freezing and thawing. Store at -20°C. Aliquots are stable for 6 months.
<b>Usage Summary:</b>	The applications listed above are for the unconjugated form of this antibody. The conjugated antibody has not been tested.
<b>Uses:</b>	IHC, Immunofluorescence, Western blot, Flow Cytometry, ELISA (Optimal dilution to be determined by the researcher)
<b>Size:</b>	200 µl
<b>Requested From:</b>	Japan

Laboratory Reagent For In Vitro Research Use Only

Not for resale without prior written consent from LifeSpan BioSciences, Inc.

Created on 9/27/2014

© 2014 LifeSpan BioSciences