

CDKN1A / WAF1 / p21 Rabbit anti-Human Polyclonal (aa116-146) (FITC) Antibody - LS-C208353 - LSBio	
CatalogID:	LS-C208353
Target:	cyclin-dependent kinase inhibitor 1A (p21, Cip1) (CDKN1A)
Synonyms:	CDKN1A Antibody, CDKN1 Antibody, CIP1 Antibody, CDK-interacting protein 1 Antibody, DNA synthesis inhibitor Antibody, p21CIP1 Antibody, MDA6 Antibody, PIC1 Antibody, SDI1 Antibody, WAF1 Antibody, CAP20 Antibody, CDK-interaction protein 1 Antibody, MDA-6 Antibody, p21 Antibody
Host	CDKN1A antibody was produced in Rabbit
Clonality:	Polyclonal
Isotype:	IgG
Conjugations:	Fluorescein (FITC)
Immunogen Species:	CDKN1A / WAF1 / p21 antibody was raised against Human
Antigen Type:	Synthetic peptide
Immunogen:	CDKN1A / WAF1 / p21 antibody was raised against kLH-conjugated synthetic peptide (aa116-146) from the C-terminus of human CDKN1A.
Specificity:	Human CDKN1A / WAF1 / p21
Epitope:	aa116-146
Reactivity:	Human, Mouse
Purification:	Protein A purified
Presentation:	PBS, pH 7.2, 0.09% sodium azide.
Recommended Storage:	May be stored at 4°C for short-term only. Aliquot to avoid repeated freezing and thawing. Store at -20°C. Aliquots are stable for 6 months.
Usage Summary:	The applications listed above are for the unconjugated form of this antibody. The conjugated antibody has not been tested.
Uses:	IHC, Western blot, ELISA (Optimal dilution to be determined by the researcher)
Size:	200 µl
Requested From:	Japan
Laboratory Reagent For In Vitro Research Use Only	
Not for resale without prior written consent from LifeSpan BioSciences, Inc.	
Created on 9/27/2014	
© 2014 LifeSpan BioSciences	