

CDKN2A / p16INK4a Rabbit anti-Human Polyclonal (aa1-30) (FITC) Antibody - LS-C208383 - LSBio	
<b>CatalogID:</b>	LS-C208383
<b>Target:</b>	cyclin-dependent kinase inhibitor 2A (CDKN2A)
<b>Synonyms:</b>	CDKN2A Antibody, ARF Antibody, Cdk4 inhibitor Antibody, CDK4 inhibitor p16-INK4 Antibody, INK4 Antibody, INK4A Antibody, Multiple tumor suppressor 1 Antibody, MTS-1 Antibody, MTS1 Antibody, p16-INK4 Antibody, p16-INK4A Antibody, p16INK4 Antibody, p19 Antibody, p14 Antibody, p16 Antibody, p16 INK4A Antibody, p16INK4A Antibody, p14ARF Antibody, TP16 Antibody, CDK4I Antibody, CDKN2 Antibody, CMM2 Antibody, MLM Antibody, p19ARF Antibody
<b>Host</b>	CDKN2A antibody was produced in Rabbit
<b>Clonality:</b>	Polyclonal
<b>Isotype:</b>	IgG
<b>Conjugations:</b>	Fluorescein (FITC)
<b>Immunogen Species:</b>	CDKN2A / p16INK4a antibody was raised against Human
<b>Antigen Type:</b>	Synthetic peptide
<b>Immunogen:</b>	CDKN2A / p16INK4a antibody was raised against kLH-conjugated synthetic peptide (aa1-30) from the N-terminus of human CDKN2A.
<b>Specificity:</b>	Human CDKN2A / p16INK4a
<b>Epitope:</b>	aa1-30
<b>Reactivity:</b>	Human
<b>Purification:</b>	Protein A purified
<b>Presentation:</b>	PBS, pH 7.2, 0.09% sodium azide.
<b>Recommended Storage:</b>	May be stored at 4°C for short-term only. Aliquot to avoid repeated freezing and thawing. Store at -20°C. Aliquots are stable for 6 months.
<b>Usage Summary:</b>	The applications listed above are for the unconjugated form of this antibody. The conjugated antibody has not been tested.
<b>Uses:</b>	Western blot, ELISA (Optimal dilution to be determined by the researcher)
<b>Size:</b>	200 µl
<b>Requested From:</b>	Japan
Laboratory Reagent For In Vitro Research Use Only	
Not for resale without prior written consent from LifeSpan BioSciences, Inc.	
Created on 9/27/2014	
© 2014 LifeSpan BioSciences	