

CLTB Rabbit anti-Human Polyclonal (Biotin) Antibody - LS-C208471 - LSBio	
CatalogID:	LS-C208471
Target:	clathrin, light chain B (CLTB)
Synonyms:	CLTB Antibody, Clathrin, light chain B Antibody, Clathrin light chain B Antibody, LCB Antibody, Clathrin, light chain (Lcb) Antibody
Host	CLTB antibody was produced in Rabbit
Clonality:	Polyclonal
Isotype:	IgG
Conjugations:	Biotin
Immunogen Species:	CLTB antibody was raised against Human
Antigen Type:	Recombinant protein
Immunogen:	CLTB antibody was raised against recombinant human CLTB.
Specificity:	Human CLTB
Reactivity:	Human
Purification:	Protein A purified
Presentation:	PBS, pH 7.2, 0.09% sodium azide.
Recommended Storage:	May be stored at 4°C for short-term only. Aliquot to avoid repeated freezing and thawing. Store at -20°C. Aliquots are stable for 6 months.
Usage Summary:	The applications listed above are for the unconjugated form of this antibody. The conjugated antibody has not been tested.
Uses:	Western blot, ELISA (Optimal dilution to be determined by the researcher)
Size:	200 µl
Requested From:	Japan
Laboratory Reagent For In Vitro Research Use Only	
Not for resale without prior written consent from LifeSpan BioSciences, Inc.	
Created on 9/27/2014	
© 2014 LifeSpan BioSciences	