

CLTB Rabbit anti-Human Polyclonal (Azide-free) (HRP) Antibody - LS-C208473 - LSBio

CatalogID:	LS-C208473
Target:	clathrin, light chain B (CLTB)
Synonyms:	CLTB Antibody, Clathrin, light chain B Antibody, Clathrin light chain B Antibody, LCB Antibody, Clathrin, light chain (Lcb) Antibody
Host	CLTB antibody was produced in Rabbit
Clonality:	Polyclonal
Isotype:	IgG
Conjugations:	Horseradish Peroxidase (HRP)
Modifications:	Azide-free
Immunogen Species:	CLTB antibody was raised against Human
Antigen Type:	Recombinant protein
Immunogen:	CLTB antibody was raised against recombinant human CLTB.
Specificity:	Human CLTB
Reactivity:	Human
Purification:	Protein A purified
Presentation:	PBS, pH 7.2. No preservative added (azide free)
Recommended Storage:	May be stored at 4°C for short-term only. Aliquot to avoid repeated freezing and thawing. Store at -20°C. Aliquots are stable for 6 months.
Usage Summary:	The applications listed above are for the unconjugated form of this antibody. The conjugated antibody has not been tested.
Uses:	Western blot, ELISA (Optimal dilution to be determined by the researcher)
Size:	200 µl
Requested From:	Japan

Laboratory Reagent For In Vitro Research Use Only

Not for resale without prior written consent from LifeSpan BioSciences, Inc.

Created on 9/27/2014

© 2014 LifeSpan BioSciences