

CLU / Clusterin Rabbit anti-Human Polyclonal (aa71-99) (PE) Antibody - LS-C208481 - LSBio	
CatalogID:	LS-C208481
Target:	clusterin
Synonyms:	CLU Antibody, APO-J Antibody, Apolipoprotein J Antibody, CLU1 Antibody, Clusterin Antibody, Complement lysis inhibitor Antibody, CLI Antibody, AAG4 Antibody, Aging-associated protein 4 Antibody, Ku70-binding protein 1 Antibody, KUB1 Antibody, SGP2 Antibody, APOJ Antibody, SGP-2 Antibody, Sulfated glycoprotein 2 Antibody, CLU2 Antibody, Complement cytolysis inhibitor Antibody, NA1/NA2 Antibody, SP-40 Antibody
Host	CLU antibody was produced in Rabbit
Clonality:	Polyclonal
Isotype:	IgG
Conjugations:	Phycoerythrin (PE)
Immunogen Species:	CLU / Clusterin antibody was raised against Human
Antigen Type:	Synthetic peptide
Immunogen:	CLU / Clusterin antibody was raised against kLH-conjugated synthetic peptide (aa71-99) from the N-terminus of human CLU.
Specificity:	Human CLU / Clusterin
Epitope:	aa71-99
Reactivity:	Human
Purification:	Protein A purified
Presentation:	PBS, pH 7.2, 0.09% sodium azide.
Recommended Storage:	May be stored at 4°C for short-term only. Aliquot to avoid repeated freezing and thawing. Store at -20°C. Aliquots are stable for 6 months.
Usage Summary:	The applications listed above are for the unconjugated form of this antibody. The conjugated antibody has not been tested.
Uses:	IHC, Western blot, Flow Cytometry, ELISA (Optimal dilution to be determined by the researcher)
Size:	200 µl
Requested From:	Japan
Laboratory Reagent For In Vitro Research Use Only	
Not for resale without prior written consent from LifeSpan BioSciences, Inc.	
Created on 9/27/2014 © 2014 LifeSpan BioSciences	