

CRKL Rabbit anti-Mouse Polyclonal (aa275-303) (FITC) Antibody - LS-C208659 - LSBio

CatalogID:	LS-C208659
Target:	v-crk avian sarcoma virus CT10 oncogene homolog-like (CRKL)
Synonyms:	CRKL Antibody, Crk-like protein Antibody
Host	CRKL antibody was produced in Rabbit
Clonality:	Polyclonal
Isotype:	IgG
Conjugations:	Fluorescein (FITC)
Immunogen Species:	CRKL antibody was raised against Mouse
Antigen Type:	Synthetic peptide
Immunogen:	CRKL antibody was raised against kLH-conjugated synthetic peptide (aa275-303) from the C-terminus of mouse Crkl.
Specificity:	Mouse CRKL
Epitope:	aa275-303
Reactivity:	Mouse
Purification:	Protein A purified
Presentation:	PBS, pH 7.2, 0.09% sodium azide.
Recommended Storage:	May be stored at 4°C for short-term only. Aliquot to avoid repeated freezing and thawing. Store at -20°C. Aliquots are stable for 6 months.
Usage Summary:	The applications listed above are for the unconjugated form of this antibody. The conjugated antibody has not been tested.
Uses:	Western blot, ELISA (Optimal dilution to be determined by the researcher)
Size:	200 µl
Requested From:	Japan

Laboratory Reagent For In Vitro Research Use Only

Not for resale without prior written consent from LifeSpan BioSciences, Inc.

Created on 9/27/2014

© 2014 LifeSpan BioSciences