

DEDD Rabbit anti-Human Polyclonal (aa114-141) (FITC) Antibody - LS-C208786 - LSBio

CatalogID:	LS-C208786
Target:	death effector domain containing (DEDD)
Synonyms:	DEDD Antibody, CASP8IP1 Antibody, DEDPro1 Antibody, DEDD1 Antibody, DEFT Antibody, FLDED-1 Antibody, FLDED1 Antibody, KE05 Antibody
Host	DEDD antibody was produced in Rabbit
Clonality:	Polyclonal
Isotype:	IgG
Conjugations:	Fluorescein (FITC)
Immunogen Species:	DEDD antibody was raised against Human
Antigen Type:	Synthetic peptide
Immunogen:	DEDD antibody was raised against KLH-conjugated synthetic peptide (aa114-141) from human DEDD.
Specificity:	Human DEDD
Epitope:	aa114-141
Reactivity:	Human
Purification:	Protein A purified
Presentation:	PBS, pH 7.2, 0.09% sodium azide.
Recommended Storage:	May be stored at 4°C for short-term only. Aliquot to avoid repeated freezing and thawing. Store at -20°C. Aliquots are stable for 6 months.
Usage Summary:	The applications listed above are for the unconjugated form of this antibody. The conjugated antibody has not been tested.
Uses:	Western blot, ELISA (Optimal dilution to be determined by the researcher)
Size:	200 µl
Requested From:	Japan

Laboratory Reagent For In Vitro Research Use Only

Not for resale without prior written consent from LifeSpan BioSciences, Inc.

Created on 9/27/2014

© 2014 LifeSpan BioSciences