

4E-BP2 / EIF4EBP2 Rabbit anti-Human Polyclonal (aa87-115) (FITC) Antibody - LS-C208870 - LSBio

CatalogID:	LS-C208870
Target:	eukaryotic translation initiation factor 4E binding protein 2 (EIF4EBP2)
Synonyms:	EIF4EBP2 Antibody, 4E-BP2 Antibody, 4EBP2 Antibody, EIF4E-binding protein 2 Antibody, PHASII Antibody
Host	EIF4EBP2 antibody was produced in Rabbit
Clonality:	Polyclonal
Isotype:	IgG
Conjugations:	Fluorescein (FITC)
Immunogen Species:	4E-BP2 / EIF4EBP2 antibody was raised against Human
Antigen Type:	Synthetic peptide
Immunogen:	4E-BP2 / EIF4EBP2 antibody was raised against kLH-conjugated synthetic peptide (aa87-115) from the C-terminus of human EIF4EBP2.
Specificity:	Human 4E-BP2 / EIF4EBP2
Epitope:	aa87-115
Reactivity:	Human
Purification:	Protein A purified
Presentation:	PBS, pH 7.2, 0.09% sodium azide.
Recommended Storage:	May be stored at 4°C for short-term only. Aliquot to avoid repeated freezing and thawing. Store at -20°C. Aliquots are stable for 6 months.
Usage Summary:	The applications listed above are for the unconjugated form of this antibody. The conjugated antibody has not been tested.
Uses:	Western blot, ELISA (Optimal dilution to be determined by the researcher)
Size:	200 µl
Requested From:	Japan

Laboratory Reagent For In Vitro Research Use Only

Not for resale without prior written consent from LifeSpan BioSciences, Inc.

Created on 9/27/2014

© 2014 LifeSpan BioSciences