

Ornithine Carbamoyltransferase Chain F Rabbit anti-E. coli Polyclonal (FITC) Antibody - LS-C210673 - LSBio

CatalogID:	LS-C210673
Target:	Ornithine Carbamoyltransferase Chain F
Host	Rabbit
Clonality:	Polyclonal
Isotype:	IgG
Conjugations:	Fluorescein (FITC)
Immunogen Species:	E. coli
Antigen Type:	Recombinant protein
Immunogen:	Recombinant E.coli Ornithine carbamoyltransferase chain F protein.
Specificity:	E. coli Ornithine Carbamoyltransferase Chain F
Reactivity:	E. coli
Purification:	Caprylic acid and ammonium sulfate precipitation
Presentation:	PBS, pH 7.4, 0.03% Proclin 300, 50% glycerol.
Recommended Storage:	Aliquot and store at -20°C or -80°C. Avoid freeze-thaw cycles.
Uses:	IHC, Western blot, ELISA (Optimal dilution to be determined by the researcher)
Size:	100 µg
Requested From:	Japan

Laboratory Reagent For In Vitro Research Use Only

Not for resale without prior written consent from LifeSpan BioSciences, Inc.

Created on 9/27/2014

© 2014 LifeSpan BioSciences