

Cystathionine Beta-lyase metC Rabbit anti-E. coli Polyclonal (FITC) Antibody - LS-C210676 - LSBio	
<b>CatalogID:</b>	LS-C210676
<b>Target:</b>	Cystathionine Beta-lyase metC
<b>Host</b>	Rabbit
<b>Clonality:</b>	Polyclonal
<b>Isotype:</b>	IgG
<b>Conjugations:</b>	Fluorescein (FITC)
<b>Immunogen Species:</b>	E. coli
<b>Antigen Type:</b>	Recombinant protein
<b>Immunogen:</b>	Recombinant Escherichia coli Cystathionine beta-lyase metC protein.
<b>Specificity:</b>	E. coli Cystathionine Beta-lyase metC
<b>Reactivity:</b>	E. coli
<b>Purification:</b>	Caprylic acid and ammonium sulfate precipitation
<b>Presentation:</b>	PBS, pH 7.4, 0.03% Proclin 300, 50% glycerol.
<b>Recommended Storage:</b>	Aliquot and store at -20°C or -80°C. Avoid freeze-thaw cycles.
<b>Uses:</b>	IHC, Western blot, ELISA (Optimal dilution to be determined by the researcher)
<b>Size:</b>	100 µg
<b>Requested From:</b>	Japan
<b>Laboratory Reagent For In Vitro Research Use Only</b> <b>Not for resale without prior written consent from LifeSpan BioSciences, Inc.</b> Created on 9/27/2014 © 2014 LifeSpan BioSciences	