

LGB / Beta-Lactoglobulin Rabbit anti-Bovine Polyclonal (FITC) Antibody - LS-C210696 - LSBio

CatalogID:	LS-C210696
Target:	progestagen-associated endometrial protein (PAEP)
Synonyms:	PAEP Antibody, LGB Antibody, Beta-Lactoglobulin Antibody
Host	PAEP antibody was produced in Rabbit
Clonality:	Polyclonal
Isotype:	IgG
Conjugations:	Fluorescein (FITC)
Immunogen Species:	LGB / Beta-Lactoglobulin antibody was raised against Bovine
Antigen Type:	Purified protein
Immunogen:	LGB / Beta-Lactoglobulin antibody was raised against bovine milk beta - Lactoglobulin.
Specificity:	Bovine LGB / Beta-Lactoglobulin
Reactivity:	Bovine
Purification:	Caprylic acid and ammonium sulfate precipitation
Presentation:	PBS, pH 7.4, 0.03% Proclin 300, 50% glycerol.
Recommended Storage:	Aliquot and store at -20°C or -80°C. Avoid freeze-thaw cycles.
Uses:	IHC, Western blot, ELISA (Optimal dilution to be determined by the researcher)
Size:	100 µg
Requested From:	Japan

Laboratory Reagent For In Vitro Research Use Only

Not for resale without prior written consent from LifeSpan BioSciences, Inc.

Created on 9/27/2014

© 2014 LifeSpan BioSciences