

TNF Alpha Rabbit anti-Bovine Polyclonal (FITC) Antibody - LS-C210700 - LSBio

CatalogID:	LS-C210700
Target:	tumor necrosis factor (TNF)
Synonyms:	TNF Antibody, APC1 protein Antibody, Cachectin Antibody, DIF Antibody, TNF Alpha Antibody, Tumor necrosis factor-alpha Antibody, TNF-alpha Antibody, TNF-a Antibody, TNFA Antibody, Tumor necrosis factor Antibody, TNF, macrophage-derived Antibody, TNF, monocyte-derived Antibody, TNFSF2 Antibody
Family / Subfamily:	TNF / not assigned-TNF
Host	TNF antibody was produced in Rabbit
Clonality:	Polyclonal
Isotype:	IgG
Conjugations:	Fluorescein (FITC)
Immunogen Species:	TNF Alpha antibody was raised against Bovine
Antigen Type:	Recombinant protein
Immunogen:	TNF Alpha antibody was raised against recombinant Bovine Tumor necrosis factor protein.
Specificity:	Bovine TNF Alpha
Reactivity:	Bovine
Purification:	Caprylic acid and ammonium sulfate precipitation
Presentation:	PBS, pH 7.4, 0.03% Proclin 300, 50% glycerol.
Recommended Storage:	Aliquot and store at -20°C or -80°C. Avoid freeze-thaw cycles.
Uses:	IHC, Western blot, ELISA (Optimal dilution to be determined by the researcher)
Size:	100 µg
Requested From:	Japan

Laboratory Reagent For In Vitro Research Use Only

Not for resale without prior written consent from LifeSpan BioSciences, Inc.

Created on 9/27/2014

© 2014 LifeSpan BioSciences