

**ALB / Serum Albumin Rabbit anti-Guinea pig Polyclonal (FITC) Antibody - LS-C210752 - LSBio**

<b>CatalogID:</b>	LS-C210752
<b>Target:</b>	albumin
<b>Synonyms:</b>	ALB Antibody, Albumin (32 AA) Antibody, Albumin (AA 34) Antibody, Albumin Antibody, Growth-inhibiting protein 20 Antibody, PRO1341 Antibody, PRO0883 Antibody, Serum albumin Antibody, PRO0903 Antibody
<b>Host</b>	ALB antibody was produced in Rabbit
<b>Clonality:</b>	Polyclonal
<b>Isotype:</b>	IgG
<b>Conjugations:</b>	Fluorescein (FITC)
<b>Immunogen Species:</b>	ALB / Serum Albumin antibody was raised against Guinea pig
<b>Antigen Type:</b>	Purified protein
<b>Immunogen:</b>	ALB / Serum Albumin antibody was raised against serum Albumin purified from guinea pig serum.
<b>Specificity:</b>	Guinea pig ALB / Serum Albumin
<b>Reactivity:</b>	Guinea pig
<b>Purification:</b>	Caprylic acid and ammonium sulfate precipitation
<b>Presentation:</b>	PBS.
<b>Recommended Storage:</b>	Aliquot and store at -20°C or -80°C. Avoid freeze-thaw cycles.
<b>Uses:</b>	IHC, Western blot, ELISA (Optimal dilution to be determined by the researcher)
<b>Size:</b>	100 µg
<b>Requested From:</b>	Japan

Laboratory Reagent For In Vitro Research Use Only

Not for resale without prior written consent from LifeSpan BioSciences, Inc.

Created on 9/27/2014

© 2014 LifeSpan BioSciences